



Be Right™



Bottle, m-endo broth, glass, 50 tests, 100 mL

Продукт №:

2373542

Опасен

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.



m-Endo broth detects total coliform bacteria (colonies with a greenish-gold sheen). For potable water samples, confirm with Lauryl Tryptose and Brilliant Green Bile Broth tubes. Samples that test positive for total coliform bacteria should be further analyzed for the presence of fecal coliform or E. coli using EC Medium or EC Medium with MUG respectively.

Спецификации

Media Type:	m-Endo Broth
Материал:	Glass
Обем / опаковка, размер:	1 each