



## Test kit, nitrogen, ammonia, model NI-SA

Продукт №:

2428700

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

**Designed for testing ammonia in freshwater or seawater**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Спецификации

|                           |                                                             |
|---------------------------|-------------------------------------------------------------|
| Smallest increment:       | 0.1                                                         |
| Бр. тестове:              | 100                                                         |
| Вид на корпуса:           | D                                                           |
| Измерване, диапазон:      | 0 - 2,4 mg/L NH <sub>3</sub> -N                             |
| Метод, име:               | Color disc/Salicylate                                       |
| Модел на тестов комплект: | NI - SA                                                     |
| Параметър:                | Nitrogen, Ammonia, mid range - For freshwater or seawater** |
| Платформа:                | Цветен диск                                                 |