



Test kit, saltwater aquaculture, model FF-3

Продукт №:

243003

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Maintaining good water quality is essential for profitable fish production. In many cases, the availability and cost of "crop" insurance is directly related to testing programs.

Designed to meet the demands of specialized salt-water applications

Each kit includes the apparatus and reagents necessary to measure 10 critical parameters using digital titrator, titration, and color disc methods

Two tests, ammonia-nitrogen and a convenient chloride titration test for salinity, have been designed exclusively for this kit

Спецификации

Модел:	FF-3
Платформа:	Цифров титратор
Приложение:	Aquaculture