



Reagent set, hardness (total), HA-71A test kit

Продукт №:

2437800

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



For determination of Total Hardness by the Drop Count EDTA Titration Method. Replacement reagent set for HA-71A Test Kit. Set contains Buffer Solution (Hardness 1), Hardness 2 and Hardness 3 Test Solutions for approximately 100 tests.

Спецификации

Измерване, диапазон:

1 - 20 mg/L CaCO₃

Метод, име:

EDTA Titration

Параметър:

Hardness (Total)

Специална функция:

Reagent set, HA-71A test kit