



Reagent set, chlorine (total), CN-65 test kit

Продукт №:

2438600

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



For determination of Total Chlorine by Drop Count Titration. Replacement reagent set for CN-65 Test Kit (225401). Set includes: Sulfamic Acid Powder Pillows, Sulfite 1 Reagent Powder Pillows and Sodium Thiosulfate Standard Solution for approximately 100 tests.

Спецификации

Измерване, диапазон:

0.2 - 4 mg/L, 1 - 20 mg/L

Метод, име:

Iodometric Titration

Специална функция:

CN-65 Test Kit