



Flask, Erlenmeyer, Wide Mouth, 2000 mL, Subdivisions 600 mL

Продукт №:

2489454

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included. Requires Size 13 rubber/silicone stopper.

Спецификации

Единица:	1 each
Капацитет:	2000 mL
Материал:	Glass Wide Mouth