



Bottle, sampling, presterilized, round, 12/pk

Продукт №:

2495012

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Works well with P-A Broth Ampules and PathoScreen Medium.

Спецификации

Dechlorinating Agent?:

No

Обем / опаковка, размер:

12/pk