



## Sample Cell: 1" Square Glass 10 mL matched set of 8

Продукт №:	2495408
BGN Цена:	Свържете се с нас
Налично	

Comes as a matched set of 8 cells.

### Use in Combination with Hach Spectrophotometers

Hach spectrophotometers, chemistries, and methods are designed, developed, and manufactured to optimize your results.

---

### Спецификации

Бележка под линия:	1 pair of cells is included with purchase of DR 5000 UV-Vis, DR 3800, DR 2800, or DR 2700 Spectrophotometer
Обем / опаковка, размер:	8 /pk