



pH 7.0 Buffer, Standard Cell Solution, 500mL for Process pH Probes

Продукт №:

25M1A1001-115 **Опасен**

BGN Цена:

Свържете се с нас

Налично



Standard Cell Solution for Hach pH/D Differential pH and ORP Sensors, packaged in resealable 500 mL bottle. Used to replenish standard cell chamber while replacing salt bridge.

Спецификации

Инструмент:	pH/D
Концентрация:	pH 7.0
Обем / опаковка, размер:	500 mL
Параметър:	pH