



Syringe filter, non-sterile PDVF, 33 mm, 0.45 μ m, 50/pk

Продукт №:	2513603
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Hach syringe filters offer an efficient way to remove particulates from samples prior to analysis. Filters feature 20% more surface area than standard 25 mm filters, which increases flow rates and capacity. PVDF (polyvinylidene fluoride) filters are ideal for solutions containing proteins. Standard Female Luer-Lok inlet with standard male slip Luer outlet. Hydrophilic. Diameter: 33mm. Pore Size: 0.45 μ m. 50/PK.

Спецификации

Диаметър:	33 mm
Единица:	50/pk
Материал:	PVDF
Обем / опаковка, размер:	50 /pk
Пори, размер:	0.45 μ m