



Parafilm M Laboratory Film

Продукт №:	251764
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Seals labware quickly, effectively. Film prevents sample evaporation and contamination; protects anhydrous materials. Can be folded repeatedly and twisted around sharp surfaces without rupturing. Stretched lengthwise at 21 °C (71 °F), Parafilm M increases approximately 200%. Becomes more flexible, softer, stronger as temperature increases. Resistant to alcohol, medium concentration acid, and alkalies. Width: 10 cm (4"); 38 m (125') roll.

Спецификации

Обем / опаковка, размер:	1 each
--------------------------	--------