



Bottle rack, plastic-coated, 15 compartments

Продукт №: 2551800

ОСТАРЯЛ АРТИКУЛ

Този артикул вече не е наличен.

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Спецификации

Обем / опаковка, размер:	1
Специална функция:	Rack, 15 Bottle Compartments