



Sample Container, 120 mL High PF

Продукт №:	2552200
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

These 120 mL (4oz.) High Profile Corning Snap-Seal™ polypropylene containers are designed for a wide variety of applications and provide a reliable leak-tight seal when closed properly. Approximately 3.5 inches tall (91 mm) and 1.75 inches (45 mm) diameter. Plastic, snap seal, translucent. Sold individually.

Спецификации

Диаметър, външен:	45 mm (approx)
Капацитет:	120 mL - 1 each
Обем / опаковка, размер:	1 /pk