



## Test kit, surface water

Продукт №:

2559800

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Опасен



The ideal tool for source water monitoring and wastewater treatment plants, this kit uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, nitrate, dissolved oxygen, phosphorus, and temperature. The rugged, lightweight kit includes the convenient pH Pocket Pal(tm) Tester and all necessary apparatus and reagents for testing.

**Ideal tool for both industrial and municipal water treatment plants**

**Uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, titrate, dissolved oxygen, phosphorus, and temperature**

**Rugged, lightweight kit includes the convenient pH Pocket Pal<sup>™</sup> Tester and all necessary apparatus and reagents**

---

### Спецификации

Бр. тестове:	100
Платформа:	Цветен диск
Приложение:	Surface Water