



Filter, TSS glass fiber, pore size 1.5 μm , diameter 110 mm, 100/pk

Продукт №:

2569100

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 μm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Спецификации

Диаметър:	110 mm
Материал:	Glass Microfiber
Обем / опаковка, размер:	100 /pk
Пори, размер:	1.5 μm