



Test kit, chloride, model CDS-DT

Продукт №:

2580600

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Спецификации

Smallest increment:	10
Бр. тестове:	100
Вид на корпуса:	F
Измерване, диапазон:	10 - 10,000 mg/L Cl ⁻
Метод, име:	Digital Titrator/ Silver nitrate
Модел на тестов комплект:	CDS-DT
Параметър:	Chloride - As Cl ⁻ , 10-10,000 mg/L Cl ⁻
Платформа:	Цифров титратор

Задължителни аксесоари

- Пипета, вариращ обем, 1,0 - 5,0 mL (Item BBP065)