



Thermometer, Enviro-Safe, Non-Mercury, -20 - 110 °C

Продукт №:

2635702

BGN Цена:

Свържете се с нас

Доставка до 2 седмици

Enviro-Safe non-Mercury thermometers feature an exclusive biodegradable liquid and dye. Easy to read, aesthetically distinctive due to dark green liquid and permanently fused markings. Each thermometer, individually calibrated against NIST traceable equipment, identified by unique serial number for ISO9000 or other quality systems.

Спецификации

Деления на термометъра:	1 °C
Обем / опаковка, размер:	1 each
Температурен диапазон:	-20 - 110 °C
Тип на флуида в термометъра:	Non-Mercury