



Cylinder, Graduated, 1000 mL, Certified

Продукт №: 2636253
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Спецификации

Вид:	Две метрични стойности
Височина:	465 mm
Измерване метод:	To Deliver
Капацитет:	1000 mL
Материал на тялото (измервателния уред):	Glass
Обем / опаковка, размер:	1 each
Сертификати:	Class A