



Buret, Class A, Glass, 50 mL

Продукт №:

2636541

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

Спецификации

Капацитет:	50 mL
Материал на тялото (измервателния уред):	Class A Glass
Специална функция:	Buret