



Be Right™



Gas washing bottle, 1200 mL

Продукт №:	2662200
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis—Hach Method 10065.

Спецификации

Body Material:	Glass
Package size UOM:	1pk
Капацитет:	1200 mL
Материал:	Borosilicate Glass
Обем / опаковка, размер:	1 each