



Test kit, hypochlorite, model CN-HR

Продукт №:

2687200

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Спецификации

Smallest increment:	0.5%
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	50 - 150 g/L Cl_2 (5-15%)
Метод, име:	Drop count titration/ Thiosulfate
Модел на тестов комплект:	CN-HR
Параметър:	Hypochlorite, high range - (Bleach) As Cl_2
Платформа:	Drop Count Titration