



## APA6000 Alkalinity Reagent 2, Indicator Solution, 1 L

|                       |                   |
|-----------------------|-------------------|
| Продукт №:            | 2696653           |
| BGN Цена:             | Свържете се с нас |
| Доставка до 2 седмици |                   |

The APA6000 Alkalinity Process analyser uses five reagents each of which is supplied in a 1 L bottle designed to provide reagents for 30 days of continuous use.

**APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample.**

**APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration.**

**APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration.**

**The fifth reagent, APA Alkalinity Cleaning Solution (269703), cleans the system during the Prime and Instrument Clean functions.**

---

### Спецификации

|                        |   |
|------------------------|---|
| Измерване, диапазон:   | N/A   |
| Инструмент:            | APA6000   |
| Параметър:             | Alkalinity  |
| Условия на съхранение: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |