



## Plastic ampules R2A, pk/50

Продукт №: 2812350  
BGN Цена: Свържете се с нас  
Не е показана дата на доставка

R2A Broth Ampules; Plastic

R2A Plastic Ampules for detection of stressed heterotrophic bacteria; pk/50.

### Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Спецификации

Media Type:	R2A
Материал:	Plastic
Обем / опаковка, размер:	50 /pk