



Vacuum filter microfunnel, pk/50.

Продукт №: 2831500
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Спецификации

Диаметър:	47 mm
Материал:	Mixed-Cellulose Ester
Обем / опаковка, размер:	50
Пори, размер:	0.45 µm