



Test kit color disc chlorine dioxide, model CLD-2

Продукт №:

2842800

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Monitor ClO₂ levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO₂ concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Спецификации

Parameter:	Chlorine Dioxide
Smallest increment:	0.1
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	0 - 2,0 ClO ₂
Метод, име:	Color Disc/DPD
Модел:	CLD-2
Модел на тестов комплект:	CLD-2 ²
Параметър:	Chlorine Dioxide
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 2.0 ppm ClO ₂