



Bottle, BOD, 300 mL, 24/pk, serialized: 1-24

Продукт №: 2898700
BGN Цена: Свържете се с нас
Не е показана дата на доставка

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included Serialized bottles have silk-screened, sequential numbers (1-24) and include a bar code.

Спецификации

| | |
|--------------------------|---------------------------|
| Body Material: | Glass |
| Package size UOM: | 24/pk, Serialized: 1 - 24 |
| Капацитет: | 300 mL |
| Материал: | Borosilicate Glass |
| Обем / опаковка, размер: | 24 /pk |