



Bottle, wash, polyethylene, wide mouth, 1000 mL, 3/pk

Продукт №: 2927204
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Polyethylene bottle with polypropylene closure and stem molded in one piece for leak-proof service. Tip can be cutback to increase flow. Available in Narrow or Wide Mouth.

Спецификации

Stem:	Transparent Cap/Stem
Еденица:	3 pk
Капацитет:	1000 mL
Материал:	LDPE
Обем / опаковка, размер:	3 pk
Цвят:	Clear