



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Продукт №:

2936500

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Спецификации

Диаметър:	47 mm mm
Материал:	Mixed-Cellulose Ester
Обем / опаковка, размер:	600
Пори, размер:	0.45 µm