



Bottle BOD, 300 mL, disposable

Продукт №:	2943101
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Спецификации

Body Material:	Plastic
Капацитет:	300 mL
Материал:	PET
Обем / опаковка, размер:	1 each