



Fixed Volume Micropipette, 100 µL

Продукт №:	2946006
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Спецификации

Капацитет:	100 µL
Обем / опаковка, размер:	1 each
Точност:	0,3 %