



## Fixed volume, Micropipette, 1000 µL

Продукт №:	2946010
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

---

### Спецификации

Капацитет:	1000 µL
Обем / опаковка, размер:	1 each
Точност:	0,5 %