



Sodium/Potassium Ionic Strength Adjustor (ISA) Powder Pills, pk/100

Продукт №:

4451569

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Ionic Strength Adjustor (ISA) for Sodium/Potassium determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Sodium/Potassium level readings. Eliminate most interferences. Pack of 100 powder pillows.

Спецификации

Package settings:	Powder Pillows
Концентрация:	ISA
Обем / опаковка, размер:	100 /pk
Параметър:	Sodium