



Nitrate Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Продукт №:

4456369

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



Ionic Strength Adjustor (ISA) for Nitrate determinations by Ion Selective Electrode (ISE) method. Powder form prevent from diluting samples. Ideal for low Nitrate level readings. Eliminate most interferences. Pack of 100 powder pillows.

Спецификации

Package settings:	Powder Pillows
Бр. тестове:	100
Измерване, диапазон:	0,05 - 14.000 mg/L NO ₃ ⁻ -N
Концентрация:	ISA
Метод, име:	Ion Selective Electrode
Обем / опаковка, размер:	100 /pk
Описание:	Nitrate Ionic Strength Adjuster Powder Pillows
Параметър:	Nitrate