



MultiTest kit case (295 x 330 x 576 mm), blue polypropylene

Продукт №:	4660800
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Спецификации

Вид на корпуса:	J
Височина:	11 in
Размери (В x Ш x Д):	295 mm x 330 mm x 576 mm