



Beaker, griffin, low form, glass, 250 mL

Продукт №:

50046H

BGN Цена:

Свържете се с нас

Налично

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Спецификации

Единица:	1 each
Капацитет:	250 mL
Материал:	Glass, Griffin, Low Form
Материал на тялото (измервателния уред):	Glass