



Beaker, griffin, low form, glass, 400 mL

Продукт №:	50048
BGN Цена:	Свържете се с нас
Доставка до 2 седмици	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Спецификации

Единица:	1 each
Капацитет:	400 mL
Материал:	Glass, Griffin, Low Form
Материал на тялото (измервателния уред):	Glass