



## Beaker, Griffin, Low Form, Glass, 600 mL

Продукт №:	50052
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Спецификации

Единица:	1 each
Капацитет:	600 mL
Материал:	Glass, Griffin, Low Form
Материал на тялото (измервателния уред):	Glass