



Beaker, griffin, low form, glass, 100 mL, 12/pk

Продукт №: 50072H
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Спецификации

Единица:	12 /pk
Капацитет:	100 mL
Материал:	Glass, Griffin, Low Form
Материал на тялото (измервателния уред):	Glass