



APA6000 Alkalinity Reagents Set

Продукт №:

6001000

Опасен

BGN Цена:

Свържете се с нас

Не е показана дата на доставка



APA6000 Alkalinity Analyser Reagent Set includes a 1 L bottle of each: APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample, APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints

Спецификации

Измерване, диапазон:

0 - 250 mg/L & 0 - 500 mg/L

Инструмент:

APA6000

Параметър:

Alkalinity

Специална функция:

Reagent Set, 3 x 1L

Условия на съхранение:

Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.