



Buffer Powder Pillows, pH 8.30, 50 mL, pk/25

Продукт №:

89868

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Buffer for preparing monochloramine standards and detecting visual alkalinity endpoints. pH 8.30, Sample Size: 50 mL sample. Pack of 25 Powder Pillows.

Fresh and reliable

Unit-dose powder pillows assure buffer freshness and eliminate the possibility of contamination.

Спецификации

Концентрация:

pH 8,3

Описание:

Buffer Powder Pillows, 25/pk

Платформа:

Прахови вложки