



Recharge kit of 4 pre-filled cartridges with premounted 2956A membranes for A1100 oxygen sensors with 33051-xT cap only

Продукт №:

2956A-AT05

Опасен

BGN Цена:

Свържете се с нас

Доставка до 2 седмици



Includes o-rings and cleaning tools. Special low concentration electrolyte.

Спецификации

Инструмент:	Orbisphere A1100-S00T O2 sensor, Orbisphere GA2x00-S00T O2 sensor
Параметър:	O2
Условия на съхранение:	Store in a dry place at room temperature