



## **Titralab AS1000 Series Sample Changer, 12 positions, 50 mL & 150 mL beakers, for two probes**

Продукт №:	AS1000.99.12150
BGN Цена:	Свържете се с нас
Доставка до 1 седмица	

### **The AS1000 Sample Changer is the perfect laboratory assistant for automating analyses on sample series.**

The Titralab AS1000 series is an automatic sample changer used in analytical laboratories in connection with AT1000 Series titrators. The AS1000 is available in four configurations of different tray capacity that holds the sample beakers. The rotating module turns the tray to change the sample position while the lifting module moves sensors and tubes up and down into each sample. The AS1000 ensures automatic stirring of samples thanks to its integrated magnetic stirrer, and rinsing of probes between each sample analysed thanks to dedicated rinse positions. The AT1000 Series/AS1000 Series setup is entirely PC-controlled with the TM1000 software to automate sample analysis series. It archives results, covering a wide variety of aqueous sample applications such as pH/Conductivity/Alkalinity and Hardness in water, Chlorides content in food products, pH/Acidity in wine and juices, Free & Total SO<sub>2</sub> in wine.

#### **Run sample analyses in total autonomy, without user intervention**

The rotating module turns the tray to change the sample position while the lifting module moves sensors and tubes up and down into each sample. Single-radius sample trays for up to 30 samples.

#### **Complete automation to free up your time**

AS1000 Sample Changer multi-parameter automation solutions release operators' time from cumbersome repetitive analyses for more valuable tasks that cannot be automated.

---

### **Спецификации**

Number of sample positions:	12
Гаранция:	24 месеца
Какво има в кутията:	AS1000 Sample Changer contains: Lifting module, Rotating module with carousel tray, support base with magnetic stirrer, probe & tubes holder, RS-232 cable with USB/Serial adapter to connect AS1000 to PC, Power supply, set of magnetic stir bars, set of PP beakers & TM1000 PC-Software with ethernet cables to connect AT1000 to local network or PC.

---

### **Какво има в кутията**

AS1000 Sample Changer contains: Lifting module, Rotating module with carousel tray, support base with magnetic stirrer, probe & tubes holder, RS-232 cable with USB/Serial adapter to connect AS1000 to PC, Power supply, set of magnetic stir bars, set of PP beakers & TM1000 PC-Software with ethernet cables to connect AT1000 to local network or PC.