



Glass Fiber Filter, Diameter 90 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Продукт №: 2546140
BGN Цена: Свържете се с нас
Не е показана дата на доставка

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Спецификации

Диаметър:	90 mm
Единица:	100/pk
Материал:	Glass Microfiber
Обем / опаковка, размер:	100 /pk
Пори, размер:	1.5 μ m