



Radiometer Analytical XR110 Reference Electrode (Red Rod, ground joint, screw cap)

Продукт №:

B20B112

BGN Цена:

Свържете се с нас

Не е показана дата на доставка

The XR110 Reference Electrode has a glass body, Red Rod Ag/AgCl reference technology system, refillable electrolyte (KCl 3M), ground joint, and porous pin junction. The XR110 reference probe is recommended for Voltammetry applications and Educational purposes. The XR110 is a screw cap electrode and requires an external connection cable for connection to a meter. (Replace Obsolete #B20B110)

Спецификации

Гаранция:	12 месеца
Какво има в кутията:	Radiometer Analytical XR110 reference electrode, 30mL KCl 3M refill solution and Operating Instruction.
Материал:	Glass
Размери (диаметър x дължина):	8 mm x 120 mm
Разтвор за пълнене:	KCl 3M
Свързка:	Screw cap S7
Съединение на сондата:	Porous plug
Температурен диапазон:	-10 to 100°C
Тип на сондата:	Laboratory reference electrode

Какво има в кутията

Radiometer Analytical XR110 reference electrode, 30mL KCl 3M refill solution and Operating Instruction.

Задължителни аксесоари

- CL111 Electrode cable FX/FIL/1m/B (Item A94L111)
- CL114 cable, 1m, for screw cap electrode FX/S7/coax, BNC plug (Item A94L114)