



## Buffer Solution (TISAB) for EZ3007 Fluoride Analyser, 2 L

Продукт №:	APPC3007-01
BGN Цена:	Свържете се с нас
Не е показана дата на доставка	

### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Спецификации

Обем / опаковка, размер:	2 L
Параметър:	Fluoride