



## EZ3004 Chloride Analyser

Продукт №: EZ3004.XXXXXXXXXX  
BGN Цена: Свържете се с нас

### Online ion-selective Chloride measurements in industrial and environmental applications

#### ISE technology for optimal analytical performance

With limited maintenance requirements and reduced reagent consumption, the EZ3000 Series are the ideal choice for a wide range of water monitoring applications where ion-selective electrodes are the preferred analytical technique. Outstanding precision and stability is guaranteed by the temperature-controlled measurement.

#### Direct, discontinuous ISE method

Contrary to separate electrodes or other analysers in the market, the EZ3000 Series does not run a continuous measurement. The principle of discontinuous ISE analysis not only enhances control over conversion of ion activity to electric potential, it also eliminates risk of cross-contamination between cycles and reduces overall consumption of reagents.

The EZ3000 Series combine unique ISE technology with a set of analysis, control and communication features in an industrial analyser mainframe with designed for the highest performance:

- Automatic direct ion-selective measurements
- Smart automatic features
- Control and communication via industrial panel PC
- Standard 4 - 20 mA signal output with alarm processing
- Communication ports supporting connectivity to Modbus
- Multiple stream analysis
- Reduced reagent consumption

There are many additional options available. Please contact Hach for more details.

---

### Спецификации

|                      |                                                                                                                |
|----------------------|----------------------------------------------------------------------------------------------------------------|
| Demineralised water: | For rinsing                                                                                                    |
| Digital outputs:     | Optional: Modbus (TCP/IP, RS485)                                                                               |
| Earth connection:    | Dry and clean earth pole with low impedance ( $< 1 \text{ Ohm}$ ) using an earth cable of $> 2.5 \text{ mm}^2$ |

|                          |                                                                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Instrument air:          | Dry and oil free according to ISA-S7.0.01-1996 quality standard for instrument air                                                                                                                |
| Precision:               | Better than 2% full scale range for standard test solutions                                                                                                                                       |
| Reagent Requirements:    | Keep between 10 - 30 °C                                                                                                                                                                           |
| Validation:              | Automatic; frequency freely programmable                                                                                                                                                          |
| Аварийен сигнал:         | 1x malfunctioning, 4x user-configurable, max. 24 VDC/0.5 A, potential free contacts                                                                                                               |
| Автоматично почистване:  | Да                                                                                                                                                                                                |
| Аналогови изходи:        | Active 4 - 20 mA max. 500 Ohm load, standard 1, max. 8 (option)                                                                                                                                   |
| Бр. потоци проби:        | 1, 2, 4 or 8                                                                                                                                                                                      |
| Гаранция:                | 24 месеца                                                                                                                                                                                         |
| Граница на откриваемост: | 2.5 mg/L                                                                                                                                                                                          |
| Дренаж:                  | Atmospheric pressure, vented, min. Ø 32 mm                                                                                                                                                        |
| Захранване:              | 100 - 240 VAC, 50/60 Hz                                                                                                                                                                           |
|                          | Max. power consumption: 120 VA                                                                                                                                                                    |
| Измерване метод:         | Discontinuous, direct measurement by combined ion-selective electrode, conform with standard methods EPA 9212 and ASTM D512-12                                                                    |
| Измерване, диапазон:     | 10 - 100 mg/L Cl <sup>-</sup>                                                                                                                                                                     |
|                          | Optional:                                                                                                                                                                                         |
|                          | 2.5 - 25 mg/L Cl <sup>-</sup>                                                                                                                                                                     |
|                          | 5 - 50 mg/L Cl <sup>-</sup>                                                                                                                                                                       |
| Какво има в кутията:     | EZ3004 Chloride Analyser, Instruction Manual, 1 x Double Bit Door Key, 1 x Mounting Brackets, 1 x ISE Chloride Electrode and 1 x empty 2.5L Reagent Container with Fittings (for Buffer Solution) |
| Калибриране:             | Automatic, 2-point; frequency freely programmable                                                                                                                                                 |
| Клас на защита:          | Protection class: Analyser cabinet: IP44 / Panel PC: IP65                                                                                                                                         |
| Кювети за проба:         | 5 minutes                                                                                                                                                                                         |
| Материал:                | Hinged part: Thermoform ABS, door: PMMA                                                                                                                                                           |
|                          | Wall section: Galvanised steel, powder coated                                                                                                                                                     |
| Метод, име:              | Ionselective                                                                                                                                                                                      |
| Модел:                   | EZ3004                                                                                                                                                                                            |
| Околна температура:      | 10 - 30 °C ±4 °C deviation at 5 - 95% relative humidity (non-condensing)                                                                                                                          |
| Параметър:               | Chloride                                                                                                                                                                                          |
| Пауза:                   | 100 - 240 VAC, 50/60 Hz                                                                                                                                                                           |
| Проба, дебит:            | 100 - 300 mL/min                                                                                                                                                                                  |
| Проба, качество:         | Maximum particle size 100 µm, < 0.1 g/L; Turbidity < 50 NTU                                                                                                                                       |
| Проба: налягане:         | By external overflow vessel                                                                                                                                                                       |
| Размери (В x Ш x Д):     | 690 mm x 465 mm x 330 mm                                                                                                                                                                          |
| Сертификати:             | CE compliant / UL certified                                                                                                                                                                       |
| Смущения:                | Bromide, sulphide, iodide, cyanide ions may interfere. Mercury must be absent. Ammonia and thiosulphate may interfere. Fats, oil, proteins, surfactants and tar.                                  |
| Тегло:                   | 25 kg                                                                                                                                                                                             |
| Температура на проба:    | 10 - 30 °C                                                                                                                                                                                        |

## Какво има в кутията

EZ3004 Chloride Analyser, Instruction Manual, 1 x Double Bit Door Key, 1 x Mounting Brackets, 1 x ISE Chloride Electrode and 1 x empty 2.5L Reagent Container with Fittings (for Buffer Solution)