



Hardness (Calcium) Reagent Set, Buret

Продукт №:

2447000-EU

Опасен

EUR Цена:

Свържете се с нас

Доставка до 2 седмици



Reagent set for the determination of Calcium Hardness by EDTA Buret Titration Method (Hach method 8222). EPA accepted method. Range: 0 - 5000 mg/L CaCO₃. Set contains CalVer 2 Calcium Indicator Powder Pillows, Potassium Hydroxide Solution (8 N) and TitraVer (EDTA) Titrant for up to approximately 100 tests.

Спецификации

EPA compliant:	Да
Бр. тестове:	100
Измерване, диапазон:	0 - 5000 mg/L CaCO ₃
Метод, име:	EDTA Titration
Метод, номер:	8222
Параметър:	Hardness (Calcium)
Специална функция:	Reagent set, Buret

Задължителни аксесоари

- Защитни ръкавици, размер L, сини, нитрил, без пудра, 50 броя (Item SM743L)