



Be Right™



## Test kit, nitrogen, nitrate color disc, model NI-14

Продукт №:

1416100-EU

Опасен

EUR Цена:

Свържете се с нас

Не е показана дата на доставка



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO<sub>3</sub> and contains 50 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Спецификации

Parameter:	Nitrogen, Nitrate
Smallest increment:	0,05
	0,5
Бр. тестове:	50
Вид на корпуса:	B
Измерване, диапазон:	0 - 1 mg/L NO <sub>3</sub> <sup>-</sup> -N
Метод, име:	Color disc/ Cadmium reduction
Модел:	NI-14
Модел на тестов комплект:	NI-14**
Най-малки инкременти, стъпки:	0,05
Параметър:	Nitrogen, Nitrate - As NO <sub>3</sub> <sup>-</sup> -N
Платформа:	Цветен диск
Стандартен обхват на измерване:	0 - 1 mg/L NO <sub>3</sub> <sup>-</sup> -N; 0 - 10 mg/L NO <sub>3</sub> <sup>-</sup> -N