



Be Right™



## Chlorine (Free and Total) Reagent Set, DPD FEAS, Digital Titrator

Продукт №:

2445300-EU

Опасен

EUR Цена:

Свържете се с нас

Не е показана дата на доставка



For determination of Chlorine (Free and Total) by Ferrous Ethylenediammonium Sulfate Digital Titration, using DPD indicator. Range: Up to 3.00 mg/L  $\text{Cl}_2$ .

Reagent Set includes: DPD Free Chlorine Reagent Powder Pillows, DPD Total Chlorine Reagent Powder Pillows, Ferrous Ethylenediammonium Sulfate Titration Cartridge, 0.00564 N. Approximately 40-100 tests.

Boric Acid free!

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Спецификации

Бр. тестове:	40 - 100
Измерване, диапазон:	0 - 3.00 mg/L $\text{Cl}_2$
Метод, име:	DPD-FEAS

---

## Задължителни аксесоари

- Защитни ръкавици, размер L, сини, нитрил, без пудра, 50 броя (Item SM743L)