



Be Right™



## TEST KIT, HA-71A HARDNESS TOTAL

Продукт №:

145201-EU

Опасен

EUR Цена:

Свържете се с нас

Доставка до 2 седмици



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model HA-71A test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizing drop count titration, this kit provides dual ranges and contains approximately 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Спецификации

Smallest increment:	1 mg/L
	1 gpg
Бр. тестове:	100
Вид на корпуса:	D
Измерване, диапазон:	1 - 20 mg/L & 17 - 342 mg/L CaCO <sub>3</sub>
Метод, име:	Drop count titration/ EDTA-powder pillows
Модел на тестов комплект:	HA-71A
Параметър:	Hardness, total - As CaCO <sub>3</sub>
Платформа:	Drop Count Titration
Тегло:	2.5 kg